# Non-Contact Infrared Human Forehead Thermometer

## Simple to operate


The non-contact infrared thermometer is designed for measuring body temperature, specifically for use on the forehead area. Requiring no skin contact, the thermometer works by collecting heat radiation emitting from the forehead. It's a more hygienic, reliable option for fever-checking compared to intrusive oral thermometers. Receive highly accurate measurements in seconds with the single click of a button. An integrated alarm sounds when raised body temperature is detected, making it easy to inform both user and individual of a potential fever.

- +/- 0.4° F (0.2° C)
- Auto shut-off after 30 seconds to save battery with power saver
- Pre-calibrated
- Commercial grade
- 1 year warranty
- Resolution: 0.1° F (0.1° C)

### **Benefits**




FDA, CE & ISO 13485 Certified



**Temp Range:** 95° F - 109.2° F (35° C - 42.9° C)



Response Time: 1 second



**Testing Range:** 2 -3 inches (5-8 cm) away from forehead







Prod No	Description	Color	Reponse Time	Accuracy	Temperature Range	Pack
THDG986	Non-contact IR Forehead Thermometer	White	1 second	+/- 0.4°F (0.2° C)	95° F - 109.2° F (35° C - 42.9°C)	1 each (60 carton)

San Jamar 555 Koopman Lane Elkhorn, Wisconsin 53121 USA T: +1.262.723.6133 F: +1.262.226.8664 info@sanjamar.com www.sanjamar.com

### Canada

3280 Bloor Street West Suite 1140 11th Floor, Centre Tower Toronto, ON M8X 2X3 CANADA T: +1.262.723.6133 F: +1.262.226.8664 canada@sanjamar.com

Schoorstraat 26a, bus I 2220 Heist-op-den-Berg BELGIUM T: +32 15 22 81 40 F: +32 15 22 81 48 emea@sanjamar.com

Europe

#### México

Av. Paseo de la Reforma No. 350 - 10 Col. Juarez Deleg. Cuauhtemoc, CP 06600 MEXICO, D.F. T:+52 (55) 3626 0772 mexico@sanjamar.com

